



# THAMNOTOXKIT F

## TEST PROCEDURE

# 1

## PREPARATION OF STANDARD FRESHWATER

- VOLUMETRIC FLASK (1 LITER)
- 5 VIALS WITH SOLUTIONS OF  
CONCENTRATED SALTS
- DISTILLED (or deionised) WATER



## 2




POUR THE 5 VIALS  
WITH CONCENTRATED SALT SOLUTIONS  
IN  $\pm$  800 ML DISTILLED WATER,  
IN THE 1 LITER VOLUMETRIC FLASK



# 3

- FILL THE FLASK TO THE 1 LITER MARK
- AERATE FOR AT LEAST 15 MINUTES








Tube with  
*Thamnocephalus platyurus*  
cysts



*Thamnocephalus*  
*platyurus* cysts




*Thamnocephalus platyurus*  
larvae

4

### HATCHING OF THE CYSTS

CYST HATCHING SHOULD BE INITIATED 20-22 HOURS PRIOR  
TO THE START OF THE TOXICITY TEST



**5**

### **1. PREHYDRATION OF THE CYSTS**

**PREPARE 20 ML “HATCHING MEDIUM” (=DILUTED STANDARD FRESHWATER)  
BY PUTTING 2,5 ML STANDARD FRESHWATER IN A GRADUATED 25 ML CYLINDER  
AND ADDING DEIONISED WATER TO THE 20 ML MARK**



# 6


**OPEN A TUBE WITH CYSTS AND  
FILL IT WITH HATCHING MEDIUM  
(approximately 1 ml)**






# 7

- CLOSE THE TUBE WITH THE STOPPER
- SHAKE THE TUBE AT REGULAR INTERVALS DURING A 30 MINUTES PERIOD






8

## 2. TRANSFER OF THE PREHYDRATED CYSTS INTO THE HATCHING PETRI DISH

EMPTY THE CONTENTS OF THE VIAL WITH PREHYDRATED CYSTS INTO A PETRI DISH





9

- MAKE SURE THAT ALL THE CYSTS ARE TRANSFERRED BY RINSING THE TUBE WITH HATCHING MEDIUM
- ADD 10 ML HATCHING MEDIUM TO THE PETRI DISH AND SWIRL GENTLY TO DISTRIBUTE THE CYSTS EVENLY






# 10

## INCUBATION OF THE CYSTS

INCUBATE THE PETRI DISH  
FOR 20-22 HOURS AT 25 °C  
UNDER CONTINUOUS ILLUMINATION  
OF MIN. 3 000 – 4 000 LUX



11

## **PREPARATION OF THE TOXICANT DILUTIONS** (e.g. a test on an effluent)

- TAKE 5 TUBES OF 10-15 ML CONTENTS AND LABEL THEM

**C1 (100), C2 (50), C3 (25), C4 (12.5), AND C5 (6.25)**





12

- ADD 5 ML STANDARD FRESHWATER TO TUBES C2, C3, C4 AND C5



# 13

**ADD 5 ML EFFLUENT SAMPLE TO  
TUBE C1 (= 100% sample)**




# 14

- ADD 5 ML EFFLUENT TO TUBE C2
- MIX THE CONTENTS OF TUBE C2 (= 50% dilution)  
WITH THE AID OF THE PIPET







15

- TRANSFER 5 ML FROM TUBE C2 TO TUBE C3
- MIX THE CONTENTS OF TUBE C3 ( = 25% dilution ) WITH THE AID OF THE PIPET






# 16

REPEAT THE SAME PROCEDURE FOR THE NEXT DILUTIONS :

- \* 5 ML FROM TUBE C3 TO TUBE C4 (= 12,5% dilution)
- \* 5 ML FROM TUBE C4 TO TUBE C5 (= 6,25% dilution)







17

## FILLING OF THE TEST PLATE

### CONTROLS

ADD 1 ML STANDARD FRESHWATER TO EACH WELL OF COLUMN 1 (WELLS A1, B1, C1, D1)





18

### TOXICANT DILUTIONS


TRANSFER 1 ML OF TEST TUBE 5 TO EACH WELL IN COLUMN 2 (WELLS A2, B2, C2, D2)



19

REPEAT THIS PROCEDURE WITH TEST TUBES 4, 3, 2 AND 1 TO FILL THE WELLS  
OF COLUMNS 3, 4, 5 AND 6 RESPECTIVELY





**20**

**TRANSFER OF THE LARVAE FROM THE HATCHING PETRI DISH  
TO THE TEST WELLS**




**21**

- PUT THE HATCHING PETRI DISH ON THE STAGE OF THE DISSECTION MICROSCOPE
- TAKE THE MICROPIPETTE LIKE A PENCIL WITH THE INDEX FINGER AND THE THUMB TO EXERT CONTROLLED PRESSURE ON THE BULB.
- SQUEEZE THE BULB GENTLY TO PROVIDE ADEQUATE SUCTION FOR PICKING UP LARVAE.





22



**TRANSFER APPROXIMATELY 50 LARVAE FROM THE PETRI DISH TO EACH RINSING WELL IN THE FOLLOWING SEQUENCE: A1 (control), A2, A3, A4, A5 AND A6 (= increasing concentrations of toxicant)**

## 23

- PUT THE MULTIWELL PLATE ON THE STAGE OF THE DISSECTION MICROSCOPE
- TRANSFER 10 LARVAE FROM RINSING WELL A1 INTO THE 3 OTHER WELLS OF COLUMN 1 ( CUPS B1, C1 AND C1 = controls)







# 24

REPEAT THE SAME TRANSFER OF 10  
LARVAE FROM RINSING WELLS A2 TO A6  
TO THE 3 WELLS OF COLUMNS 2 TO 6  
( in this sequence, i.e. from the lowest to the  
highest toxicant concentration )





**25**

**PUT THE PARAFILM STRIP ON TOP OF THE MULTIWELL PLATE AND PUT THE COVER ON TOP**



# 26

**PUT THE MULTIWELL PLATE IN THE  
INCUBATOR AT 25 °C, IN DARKNESS,  
FOR 24 HOURS**



**- PUT THE MULTIWELL PLATE ON THE STAGE OF THE DISSECTION MICROSCOPE**

- CHECK THE WELLS OF ROWS B, C AND D AND COUNT THE NUMBER OF DEAD LARVAE IN EACH CUP
- SCORE THE MORTALITY DATA ON THE “RESULTS SHEET”
- CALCULATE THE 24h LC50 WITH AN APPROPRIATE PROGRAM