



ALGALTOXKIT M

PRODUCT LIST

Order Code

TK42 **ALGALTOXKIT M with *Phaeodactylum tricoratum*** (estuarine/marine tests)
(2 tests per Toxkit)

Including : 2 tubes with *Phaeodactylum tricoratum* algal inoculum, 2 trays with 18 long cells, 2 algal stock cells, 1 calibration long cell, 1 set of concentrated media, Standard Operating Procedure, Bench Protocol, QC Specification sheet, results sheet

MATERIALS FOR SEPARATE ORDERS

TEST BIOTA

TB42 5 tubes with *Phaeodactylum tricoratum* inoculum (estuarine/marine tests)

TEST CONTAINERS

TC41 50 long cells with lid

TC41B 270 long cells with lid (bulk)

TC41C 500 long cells with lid (bulk)

TC42 5 holding trays for long cells

TC43 5 holding trays with long cells and lids

TEST MEDIA

TM32 **5 sets of concentrated media and nutrients for marine algal tests** (for 5 x 1 l culturing medium)

TOXKIT EQUIPMENT

TE15 Spectrophotometer (Jenway 6300) for algal and protozoan tests

TE21 All-round Toxkit incubator (135 l. with 4 built-in lights)

TE21B All-round Toxkit incubator (255 l. with 6 built-in lights)

TE21C All-round Toxkit incubator (455 l. with 8 built-in lights)

TE31 Spare lamp (holder + light tube) for 135 l. & 255 l. All-round Toxkit incubator

TE32 Spare light tube for 135 l. & 255 l. All-round Toxkit Incubator

TE34 Spare light tube for 455 l. All-round Toxkit Incubator

TE35 Led lamp for 135 l., 255 l. & 455 l. All-round Toxkit incubator

TE41 Spare light bulb for Jenway 6300 Spectrophotometer