

DUCKWEED LIGHT UNIT

FOR PERFORMANCE OF THE DUCKWEED TOXKIT F




Table top Light Unit for incorporation in an incubator or temperature controlled room at 25 °C.
Uniform illumination with full spectrum grow lights for optimal growth of the duckweed.

Suited for germination of the *Spirodela polyrhiza* turions and illumination of the test plates.

Up to four complete *Spirodela* duckweed toxicity tests can be performed at once.

ORDER CODE: TE51

Specifications

Dimensions	50 cm Long x 20 cm Wide x 22 cm High
Power requirements	220-240 V AC 50 Hz
Light source	F11T5 11W HD (6400 K)



MICROBIOTESTS INC - Kleimoer 15 - 9030 Gent - BELGIUM

T: +32 9 380 85 45

info@microbiotests.com

www.microbiotests.com