


# LIGHT TABLE

FOR PERFORMANCE OF MICROBIOTESTS DAPHTOXKIT F TESTS



ORDER CODE: TE14

The light table, LightPad 930 LX light box by Artograph, is a practical tool for the performance of daphnia toxicity tests with the Daphtoxkit F.

The light table is perfectly suited for handling and scoring of the daphnias in the multiwell test plates.

## Features:

- ◆ Light box with a 30,5 × 22,9 cm lit surface
- ◆ Super bright LED lamps provides a perfectly even illuminated surface
- ◆ Variable brightness from 500 lux to super bright 5.000 lux
- ◆ Long-life LED lights last up to 50.000 hours
- ◆ LED lights color temperature of 6.500 K, cool white
- ◆ Durable and attractive extruded aluminium frame
- ◆ Double layered firm work surface
- ◆ Power supply included; AC Adapter 100 V - 240 V, 50/60 Hz



**MICROBIOTESTS INC** - Kleimoer 15 - 9030 Gent - BELGIUM

T: +32 9 380 85 45

[info@microbiotests.com](mailto:info@microbiotests.com)

[www.microbiotests.com](http://www.microbiotests.com)