

STEREO MICROSCOPE, VISISCOPE® 250

FOR PERFORMANCE OF MICROBIOTESTS TOXKIT TESTS



ORDER CODE: TE16

Reliable microscope for educational and research purposes.

SZB 250 has a 0,7 x - 4,5 x zoom range (6,43:1 zoom ratio), which is ideal for routine inspections and performing Microbiotests Toxkit tests. Supplied ready-to-use with eyepieces, objective and multiplug power supply.

- ◆ Compact and slim unit thanks to V-shaped Greenough Optical System
- ◆ Incident and transmitted light with two efficient LEDs
- ◆ Efficient LED energy saving lamp
- ◆ Adjustable angle of the incident illuminator
- ◆ Integrated handle for easy lifting and carrying of the microscope

Specifications

Model	SZB250
Observation method	Brightfield
Head	45° inclined, 360° rotating
Eyepieces	WF 10x/21 mm
Magnification	0,7x - 4,5x
Stand	Fixed, with focus and handle
Focusing system	Coaxial coarse and fine knobs
Illumination	LED (3W)
Working distance	85 mm
Interpupillary adjustment	51 - 75 mm



MICROBIOTESTS INC - Kleimoer 15 - 9030 Gent - BELGIUM

T: +32 9 380 85 45

info@microbiotests.com

www.microbiotests.com