



## RAPIDTOXKIT F THAMNO

### PRODUCT LIST

#### Order Code

TK37 **RAPIDTOXKIT F Thamno with *Thamnocephalus platyurus*** (freshwater tests)  
**(up to 45 tests per Toxkit)**

Including : 3 tubes with *Thamnocephalus platyurus* cysts, 1 flask with prepared medium, 3 hatching vessels, 6 observation plates, 6 holders with test tubes and stopper, 1 vial with red microspheres, 1 vial with fixative, Standard Operating Procedure, Bench Protocol, QC Specification sheet, results sheets

#### MATERIALS FOR SEPARATE ORDERS

##### TEST BIOTA

TB31 10 tubes with *Thamnocephalus platyurus* cysts (freshwater tests)

##### TEST MEDIA

TM41 5 flasks with prepared Standard Freshwater for Rapidtoxkit tests (5 x 150 ml EPA medium)

##### TEST ITEMS

TI14 10 hatching vessels for Rapidtoxkit tests  
TI23 5 observation plates for Rapidtoxkit tests  
TI28 5 holders with test tubes for Rapidtoxkit tests (9 tubes with stopper per holder)  
TI29 50 test tubes with stopper for Rapidtoxkit tests  
TI30 5 test tube holders for Rapidtoxkit tests  
TI31 5 vials with red microspheres for Rapidtoxkit tests  
TI31B 1 vial with red microspheres for Rapidtoxkit tests  
TI31C 3 vials with red microspheres for Rapidtoxkit tests  
TI32 5 vials with fixative for Rapidtoxkit tests

##### TOXKIT EQUIPMENT

TE16 Dissection microscope (SZB 250)  
TE21 All-round Toxkit incubator (135 l. with 4 built-in lights)  
TE21B All-round Toxkit incubator (255 l. with 6 built-in lights)  
TE21C All-round Toxkit incubator (455 l. with 8 built-in lights)  
TE22 Ring light with holder for Rapidtoxkit tests  
TE26 Ring light (without holder) for Rapidtoxkit tests  
TE27 Holder for ring light for Rapidtoxkit tests  
TE31 Spare lamp (holder + light tube) for 135 l. & 255 l. All-round Toxkit incubator

- TE32 Spare light tube for 135 l. & 255 l. All-round Toxkit Incubator
- TE34 Spare light tube for 455 l. All-round Toxkit Incubator
- TE35 Led lamp for 135 l., 255 l. & 455 l. All-round Toxkit incubator