



THAMNOTOXKIT F

PRODUCT LIST

Order Code

- TK31 **THAMNOTOXKIT F with *Thamnocephalus platyurus* (freshwater tests) (6 tests per Toxkit)**
Including : 6 tubes with *Thamnocephalus platyurus* cysts, 6 multiwell test plates, 1 sets of concentrated media,
6 micropipettes, 2 small petri dishes, 6 sheets Parafilm, Standard Operating Procedure,
Bench Protocol, QC Specification sheet, results sheets

MATERIALS FOR SEPARATE ORDERS

TEST BIOTA

- TB31 10 tubes with *Thamnocephalus platyurus* cysts (freshwater tests)

TEST CONTAINERS

- TC31 12 multiwell plates for *Thamnocephalus*, *Artemia* or *Ceriodaphnia* tests

TEST MEDIA

- TM11 10 sets of concentrated media for freshwater tests (for 10 x 1 l EPA medium)

TEST ITEMS

- TI11 10 micropipettes for rotifer and *Thamnocephalus* tests
TI21 15 small petri dishes (Ø 55 mm)
TI27 10 sheets Parafilm for rotifer and crustacean test plates (3 x 2 squares)

TOXKIT EQUIPMENT

- TE16 Dissection microscope (SZB 250)
TE21 All-round Toxkit incubator (135 l. with 4 built-in lights)
TE21B All-round Toxkit incubator (255 l. with 6 built-in lights)
TE21C All-round Toxkit incubator (455 l. with 8 built-in lights)
TE31 Spare lamp (holder + light tube) for 135 l. & 255 l. All-round Toxkit incubator
TE32 Spare light tube for 135 l. & 255 l. All-round Toxkit Incubator
TE34 Spare light tube for 455 l. All-round Toxkit Incubator
TE35 Led lamp for 135 l., 255 l. & 455 l. All-round Toxkit incubator