



ALL ROUND TOXKIT INCUBATOR

ORDER CODE TE21 / TE21B / TE21C



Thermostatically controlled incubator with built-in lights
(for all Microbiotests Toxkit tests)

- **order code TE21**
All-round Toxkit incubator (135 litres with 4 built-in lights)
Incl. 2 metal racks + 4 sockets
- **order code TE21B**
All-round Toxkit incubator (255 litres with 6 built-in lights)
Incl. 3 metal racks + 7 sockets
- **order code TE21C**
All-round Toxkit incubator (455 litres with 8 built-in lights)
Incl. 4 metal racks + 9 sockets

The All-round Toxkit incubators are used for continuous temperature control over a range of 2 °C to 40°C with or without illumination.

This makes them ideal for application as incubator for all Microbiotests Toxkits.

- ◇ **Temperature range 2 °C to 40 °C, continuously adjustable in steps of 0,1 °C**
- ◇ **Depending on the model 4, 6 or 8 built-in lights**
- ◇ **Low power consumption**
- ◇ **Illuminated display of preset/current temperature**
- ◇ **Power sockets inside the incubator**
- ◇ **3 models available: 135, 255 and 455 litres**



MICROBIOTESTS

MICROBIOTESTS INC - Kleimoer 15 - 9030 Gent- BELGIUM

T: +32 9 380 85 45

info@microbiotests.com

www.microbiotests.com

The temperature can be set in steps of 0.1 °C and an LED display shows both the set temperature and the current temperature in the cabinets.

Other electrical devices, which require a power supply, can be connected to sockets incorporated in the interior of the cabinet.

Sockets CEE 7/5, type E with hinged lid, 230 V / 16 A, 2p + E, IP 44

The integral temperature control unit meets the requirements of the EMC directive issued as IEC 61326: "Electrical devices for measurement, monitoring and for use in laboratories".

Improved, robust, insulated housing and highly efficient components provide maximum energy efficiency.

There are 3 models available with standard doors from 135 to 455 litres net capacity.

TECHNICAL DETAILS

capacity	135 Ltr. (255 or 455 Ltr.)
overall dimensions	600 deep x 600 wide x 850 high (mm) (600 x 610 x 1640 or 750 x 730 x 1640)
inside dimensions	513 deep x 441 wide x 702 high (mm) (470 x 440 x 1452 or 600 x 560 x 1452)
weight	42 kg (52 or 82 kg)
construction	fully insulated cabinet with universal plug-in temperature
control unit operation	robust, water protected foil front panel 2 keys with tactile acknowledgement
range of adjustment	2°C to 40°C; adjustable in 1°C-increments
ambient temperature	10°C - 35°C
display	3 digit for actual temperature, resolution 0,1°C 2 digit for temperature set-point
refrigerating unit	compressor-system
heating unit	heating element
fan	radial, discharge capacity 120 m3 /h
power requirement	230V/50Hz
light unit	4 (6 or 8) LED lamps (32.2 cm long) of 3.8 W that can be used individually
power supply inside	3 (or 4) extra individual plug-sockets are available



MICROBIOTESTS INC - Kleimoer 15 - 9030 Gent- BELGIUM

T: +32 9 380 85 45

info@microbiotests.com

www.microbiotests.com